


# Guide Wire



## PRODUCT FEATURES

- Nitinol core wire guarantees outstanding kink resistance (Stainless steel core wire for option)
- Teflon coating and smooth remote design with excellent integration and accessibility, provides excellent crossibility for other interventional devices
- Customized Diameters and lengths are available

## ORDERING INFORMATION

Dimensions	
OD Range	0.018"-0.038"
Length Range	45cm-300cm

0.025"	0.035"	0.038"	Length (cm)	Tip Shape
BM-GWC025-80S	BM-GWC035-80S	BM-GWC038-80S	80	Straight
BM-GWC025-80J	BM-GWC035-80J	BM-GWC038-80J	80	J
BM-GWC025-120S	BM-GWC035-120S	BM-GWC038-120S	120	Straight
BM-GWC025-120J	BM-GWC035-120J	BM-GWC038-120J	120	J
BM-GWC025-150S	BM-GWC035-150S	BM-GWC038-150S	150	Straight
BM-GWC025-150J	BM-GWC035-150J	BM-GWC038-150J	150	J
BM-GWC025-180S	BM-GWC035-180S	BM-GWC038-180S	180	Straight
BM-GWC025-180J	BM-GWC035-180J	BM-GWC038-180J	180	J
BM-GWC025-260S	BM-GWC035-260S	BM-GWC038-260S	260	Straight
BM-GWC025-260J	BM-GWC035-260J	BM-GWC038-260J	260	J