


# Radial Artery Tourniquets



## INTRODUCTION

The device is intended for use with hemostasis after the radial arterial catheterization.

## PRODUCT FEATURES

- Unique designed for radial artery hemostasis of post percutaneous puncture
- Superior fixed plate without defects and transparent material ensures clear visual control
- Soft and comfortable adjustable band
- Specifications meet the needs of clinical application

## PRODUCT COMPOSITION

The device is consist of the fixed plate, adjustable band and elastomer.

Catalogue No.	Product description
BM-RAT1	Card button
BM-RAT2	Rotary