



# Thoracic Catheter w/Trocar

Thoracic trocar catheter is a flexible plastic tube that is inserted through the side of the chest into the pleural space. It is used to remove air or fluid, or pus from the intrathoracic space.

## FEATURES

- ◆ Soft, kink-resistant smooth PVC material.
- ◆ Depth markings assist the accuracy for catheter placement.
- ◆ Smooth side eyes and open distal tip.
- ◆ Radiopaque stripe confirms placement easily.

### *Item Code*

Item Code	Size	Type
<i>Sharp Tip</i>		
I12001	8 Fr	Trocar Type
I12002	10 Fr	Trocar Type
I12003	12 Fr	Trocar Type
I12004	14 Fr	Trocar Type
I12005	16 Fr	Trocar Type
I12007	20 Fr	Trocar Type
I12008	22 Fr	Trocar Type
I12009	24 Fr	Trocar Type
I12011	28 Fr	Trocar Type
I12013	32 Fr	Trocar Type