



# Thoracic Catheter

Thoracic trocar catheter is a flexible plastic tube that is inserted through the side of the chest into the pleural space. It is used to remove air or fluid, or pus from the intrathoracic space.

## FEATURES

- ◆ Soft, kink-resistant smooth PVC material.
- ◆ Depth marking for identification of the tube position.
- ◆ Smooth side eyes and open distal tip.
- ◆ Radiopaque stripe confirms placement easily.

### *Item Code*

Item Code	Size	Tube Type
I12108	8 Fr	Straight
I12110	10 Fr	
I12112	12 Fr	
I12114	14 Fr	
I12116	16 Fr	
I12120	20 Fr	
I12122	22 Fr	
I12124	24 Fr	
I12128	28 Fr	
I12132	32 Fr	
I12034	34 Fr	
I12036	36 Fr	
I12225	16 Fr	
I12227	20 Fr	
I12229	24 Fr	
I12231	28 Fr	
I12233	32 Fr	
I12235	36 Fr	