






# Reusable Pressure Infusion Bag



## Product Advantage

Photo	Feature	Benefit
See picture above	Single patient use	Reduce the economic cost of hospitals
	Large inflatable ball, soft, good elasticity	Easier and more efficient inflation
	Unique design, equipped with Robert clip	Secondary pressure keeping to avoid air leakage, more safe and reliable
	Metal pressure indicator mark	Precise control of pressure and flow
	Transparent nylon mesh material	Not easy to cut and tear, can clearly observe the infusion bag and residual volume, convenient for rapid setting and replacement of infusion bag
	Unique hook design	To avoid the risk of falling off after the bag volume is reduced, and it is safer to use

## Scope of Application

- ★ Used to continuously pressurize the liquid containing heparin to flush the built-in arterialpiezometer.
- ★ Arterial transfusion during surgery and emergency situations.
- ★ In cerebrovascular surgery interventional treatment, in order to prevent blood from flowing back to the coaxial catheter leading to thrombosis and falling off, causing vascular embolism in the body, it is necessary to use an infusion pressurized bag for high-pressure perfusion, and continuous saline drip irrigation in the catheter.
- ★ Emergency fast blood transfusion in the field, battlefield, hospital and other occasions.

## Product Parameter

Capacity	Medlinket Ref No.	Specification	Packing	Price Code	Photo
500ml	Y000R05	Metal pressure indicator, silicone tube, material: nylon, waterproof	24pcs/box	B6	C
1000ml	Y000R10	Metal pressure indicator, silicone tube, material: nylon, waterproof	18pcs/box	B0	B
3000ml	Y000R30	Metal pressure indicator, silicone tube, material: nylon, waterproof	12pcs/box	B0	A

## 【Related Product Recommendations】 Disposable Pressure Infusion Bag



### Shenzhen Med-link Electronics Tech Co., Ltd.

Address: Zone A of 1st and 2nd Floor, and 3rd Floor, Building A, No. 7, Tongsheng Industrial Park Road, Shanghenglang Community, Dalang Street, Longhua District, 518109 Shenzhen, PEOPLE'S REPUBLIC OF CHINA  
 Tel: +86 400-058-0755 / +86 755-6112-0299  
 Email: marketing@med-linket.com Website: www.med-linket.com



\*Declaration: Shenzhen Med-Link Electronics Tech Co., Ltd. ("MedLinket") hereby declares that all registered trademarks, product names, model numbers, etc., mentioned in the document belong to respective original owners or manufacturers. The information referenced herein is intended to clearly describe the compatibility features of MedLinket products for reference only and should not be construed as an operational guideline for medical institutions or related entities. If any adverse consequences result from misunderstanding or incorrect use of the information, MedLinket shall not be held liable.