




Disposable Trocar Semi-Detachable Trocar

- One Hand Detach
- Precision
- Light Weight Easy Operation



Kit




Factory : 101-201,Bldg 38,Xialang Industrial Area,Heshuikou,Matian Street, Guangming District,Shenzhen,China

+86 755-8258 6700 +86 755-8258 6704 www.unimicromed.com marketing@unimicromed.com

Disposable Trocar Semi-Detachable Trocar



One Hand Operation
Upward Pushing Switch to Detach



Stainless Steel Obturator
Enhance Penetration of Trocar



Barb Thread Design
Innovative Design,
Stronger Fixation



Visible Tip
Unfettered Visualization During Insertion



Tip	Product Code	Specification			
Auto-Locking	UNTR-I	5*70mm			
		5*100mm	10*100mm	11*100mm	12*100mm
		5*150mm	10*150mm	11*150mm	12*150mm
Bladeless	UNTR-II	5*70mm			
		5*100mm	10*100mm	11*100mm	12*100mm
		5*150mm	10*150mm	11*150mm	12*150mm
Visible	UNTR-III	5*70mm			
		5*100mm	10*100mm	11*100mm	12*100mm
		5*150mm	10*150mm	11*150mm	12*150mm
Pediatric	UNTR-IV	3*60mm	4*60mm	5*60mm	