

Disposable Hemoclip

Indications

It is suitable for clamping the soft tissue in the digestive tract under the guidance of clinical endoscope.

Features

- ※ The front end is exquisitely designed to provide excellent passage of the clamp and smooth release under the endoscope
- ※ The clip is pre-installed on the front end of the flexible releaser, and the product consists of a single-use clip part and a flexible delivery device part
- ※ It must be used together with gastrointestinal endoscopes such as gastroscopes and colonoscopes with a biopsy forceps channel size greater than or equal to 2.8mm
- ※ It can be opened and closed repeatedly, and the function of repeated pre-clamping can be realized
- ※ Multiple span options, suitable for closure of different wounds and bleeding points
- ※ Rotatable function helps doctors grasp tissue at the most appropriate angle