

Disposable Retrieval Pouch


- Retrieval Pouch with Different Sizes
- Flexible Memory Ring and Introducer
- High Intensity TPU Pouch



Durable, Tear Resistant



Easy Insertion



Different Pouch Sizes




Disposable Retrieval Pouch



Auto Retrieval Pouch (MEBD)

The pouch can be un-cinched.
Red button design allows one-click detach.



Active Channel Retrieval Pouch (MEBC)

Be opened multiple times.
Enhanced metal ring ensures large pouch opening when loading specimen.



Retrieval Pouch with Introducer (MEBS)

Simple kind, equipped with a string.
Red bead design prevents string falling into cavity.
Equipped with a string.

Production Description	Size	200ml	400ml	600ml
Active Channel Retrieval Pouch	10mm*33cm	MEBC200	MEBC400	MEBC600
Auto Retrieval Pouch	10mm*33cm	MEBD200	MEBD400	
Retrieval Pouch with Introducer	10mm*20cm	MEBS200	MEBS400	MEBS600
SOLD 10 UNITS PER BOX				